



Starting a glioma database - The TUM experience

Benedikt Wiestler
Dept. of Neuroradiology
TU Munich



Klinikum rechts der Isar
Technische Universität München



Neuro-Kopf-Zentrum

Klinikum rechts der Isar
Technische Universität München



Why would you want a database?

- **Larger n for discovery / validation**
- **Test generalizability across scanners / centers**
- **Serve as data source for challenges**
- **...**

From Excel to a database

- **Local, with neurology / neurosurgery / neuropathology**
- **„Data information system“**
- **Some automatic information (Tumor volume, OS ...)**




...to a multi-center database

- **Bjoern Menze (TUM Computer Science)**
- **Roland Wiest (Bern Neuroradiology)**
- **Several interested centers in D-A-CH**



What data do we want?

- **Clinical?**
- **Pathology data?**
- **Imaging data / sequences?**
- **(Anonymization / Pseudonymization)**





	3D T1w Pre	Ax 2D FLAIR	Ax 2D DWI		Ax 2D T2w	3D T1w Post ^b
Sequence	IR-GRE ^{d,e}	TSE ^c	EPI ^f	Contrast Injection ^a	TSE ^c	IR-GRE ^{d,e}
Plane	Sagittal/axial	Axial	Axial		Axial	Axial/sagittal
Mode	3D	2D	2D		2D	3D
TR [ms]	2100 ^g	>6000	>5000		>2500	2100 ^g
TE [ms]	Min	100–140	Min		80–120	Min
TI [ms]	1100 ^h	2500				1100 ^h
Flip angle	10°–15°	90°/≥160°	90°/180°		90°/≥160°	10°–15°
Frequency	256	≥256	128		≥256	256
Phase	256	≥256	128		≥256	256
NEX	≥1	≥1	≥1		≥1	≥1
FOV	256 mm	240 mm	240 mm		240 mm	256 mm
Slice thickness	1 mm	3 mm	3 mm		3 mm	1 mm
Gap/apacing	0	0	0		0	0
Diffusion options			$b = 0, 500, \text{ and } 1000 \text{ s/mm}^2$ ≥3 directions			
Parallel imaging	Up to 2x	Up to 2x	Up to 2x		Up to 2x	Up to 2x
Scan time (approx)	5–8 min	4–5 min	3–5 min		3–5 min	5–8 min



Anonymization

- **Easier, „safer“**

Pseudonymization

- **Follow-up over time**

Local pre-processing software

- **DICOM -> NiFTI**
- **Co-registration, skull-stripping (saves transforms)**
- **Pseudonymization**
 - **Hash algorithm using a local password?**
 - **Local pseudonymization (list-based)?**
- **Metadata saved as .JSON**



elephant-client

Logo Bratum

1

import

2

meta data

3

3D computations

4

5

Step: 0

Scan results

We found complete nifti files for the following scans:

- Ex1
- Ex2

Shall we proceed with these scans?

► Yes, that looks good. Let's continue!

No, I need to make changes. Please take me back!

Defaced

exporter

Program Status, Backend Status



elephant-client

Logo Bratum

import

meta data

3D computations

review

export

Step: 2

Computation in Progress

Yes, that looks good. Let's continue!

exam

Ex1

Ex2

status


robex skullstripping


robex skullstripping

Show results!

No, I need to make changes. Please take me back!

Program Status, Backend Status





Central storage

- **Secure upload**
- **Central post-processing (e.g. BRATS algorithms)**
- **Encrypted storage on LRZ servers**

Data use

- **For members of the working group**
- **Central request**
- **„Opt-out“**
- **(Anonymized data for challenges)**

